Corrosion and material compatibility tests with Gigasept Instru AF use-solutions - Summary of Results -

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
Aluminium

Plastics

PE, hard	Polyethylene, hard
PS	Polystyrene
PS, glass-clear	Polystyrene, glass clear
ABS	Acrylonitrile butadiene styrene polymer
PMMA	Polymethyl methacrylate
PC	Polycarbonate
POM	Polyoxymethylene
PET	Polyethylene terephthalate
PA 6	Polyamide 6
PSU	Polysulphone
PVC, hard	Polyvinyl chloride, hard


Rubbers

NR/SBR	Natural rubber/ styrene butadiene rubber
EPDM	Ethylene propylene diene
NBR	Nitrile rubber
IIR	Isobutene-isoprene rubber
CR/SBR	Chloroprene rubber
FKM	Fluorocarbon rubber
PVC, soft	Polyvinyl chloride, soft
SI	Polydimethylsiloxane

Test method

Corrosion test

Standard test pieces (20 x 100 x 1mm) with an overall surface of 0.00424m² were 60% immersed in Gigasept Instru AF use-solutions (1.5% and 3%) at a temperature of 35°C over a 24-hour period. Any resulting weight differences (increase/decrease) were measured in g/m² and the test pieces and solutions were visually inspected.

Material compatibility

Standard test pieces (40 x 60mm) were immersed in Gigasept Instru AF use-solutions (1.5% and 3%) at a temperature of 40°C over a two-week period. Any resulting weight difference (swelling) was expressed as a % and the test materials and working solutions were visually inspected.

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Norderstedt, Germany 23.11.2004

SCHÜLKE & MAYR GMBH
Research & Development
Product Development 1

Michael Mohr