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
EDMax Absorbent Paper Points and Gutta Percha Points  
Class IIa medical devices – CE 0197

Lengths: 21-25-31 mm (except opener: 15 mm)

## Recommended Protocol

Crown Down • Speed 500 rpm • Torque Limit 1.5 N.cm • Continuous rotation • Pecking and wall brushing motion • Recapitulate • Constantly irrigate the canal and wipe the file • Stop the preparation as soon as you meet resistance at the working length

Type of canal	Opener	Glide path	Shaper				Finisher	
	20/.10	15/.03	20/.04	20/.06	25/.04	25/.06	30/.04	40/.04
<b>Large</b>	●	●	○	○	●	○	●	◐
<b>Medium</b>	●	●	○	●	○	●	○	○
<b>Narrow</b>	●	●	●	○	◐	○	○	○



For retreatments, EDMax 20/.06 also allows to efficiently remove Gutta Percha to desobturate the canal.

## The EDM Inventor for Endodontics

EDM (Electric Discharge machining) consists of high-power and high frequency electric sparks between a hair-thick metal wire cutting tool and the in-process NiTi file, that leads to locally melting and evaporating the material without mechanical contact between the electric wire and the file. This prevents the formation of microcracks in the file structure, and brings unique features such as combined core softness and surface hardness + roughness resulting in NiTi files with SOFT POWER! EDMax also offers huge potential of amazing designs such as instruments with split blades for a wide range of possible breakthrough innovative applications. Neolix, a French startup company founded in 2009, has been fully dedicated from the first day to its patented EDM process applied to Endodontics, and is the only worldwide manufacturer to master EDM for Endodontics, proudly as a pioneer.

After Neoniti file range launched in 2014 and already well established for its superior cyclic fatigue resistance, EDMax propels Neolix NiTi files up to new heights thanks to its optimized features and improved key advantages: **resistance to cyclic fatigue and torsional stress, high cutting efficiency, and respect of natural anatomy.**

**Thanks to its unique EDM process, Neolix opens a new era in Endodontics!**

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©communication - november 2021 - rev.f