

**Ricerca e Selezione
Qualità & Soluzioni
per voi!**


ilic SINCE 1949
QUALITY & VALUE

miniKUT[®]
MINIMALLY INVASIVE FILE SYSTEM




Mentre pensiamo di preservare la dentina naturale per garantire un restauro duraturo, il nome miniKUT è diventato proprio appropriato. L'obiettivo è potenziare il clinico senza limitazioni. Spero che troviate il sistema miniKUT Rotary NiTi tanto vantaggioso per i vostri pazienti quanto lo è per me...

Dr. L. Stephen Buchanan




miniKUT Shaping Files

(.03 Taper and .05 Taper - Tip sizes 20-45 - 21,25, and 31mm)




miniKUT EZP Negotiating Files

(.03 Taper and .05 Taper - Tip sizes 15 - 21,25, and 31mm)




miniKUT MB2 Files

(.05 Taper - Tip size 15 - 17mm)




miniKUT Orifice Shaping Files

(.07 Taper - 17mm - Tip size 15)



miniKUT Retreatment Files

(.06 Taper and .08 Taper - Tip size 20(06), 25(08) - 21 and 25mm)



Calibration Legend

21mm: 18, 19, 20mm

25mm: 18, 19, 20, 22, 24mm

31mm: 24, 25, 26, 28, 30mm

IVAN ILIC' srl - Via L. Muratori, 46/11 - 20135 Milano - T (+39) 02 36512990 - ilic@ilicdental.com

www.ilicdental.com



miniKUT Rotary NiTi Endodontic Shaping Files (6 Packs)

Tip Size	.03 Taper			.05 Taper		
	21mm	25mm	31mm	21mm	25mm	31mm
20	1032021	1032025	1032031	1052021	1052025	1052031
25	1032521	1032525	1032531	1052521	1052525	1052531
30	1033021	1033025	1033031	1053021	1053025	1053031
35	1033521	1033525	1033531	1053521	1053525	1053531
40	1034021	1034025	1034031	1054021	1054025	1054031
45	1034521	1034525	1034531	1054521	1054525	1054531



miniKUT Procedure Packs (3 Packs)

15/0.7, 15/0.3, 25/0.3 (17mm & 25mm)	1430325
15/0.7, 25/0.5, 35/0.5 (17mm & 25mm)	1430525



miniKUT EZP Rotary NiTi Endodontic Negotiating Files (6 Packs)

Tip Size	.03 Taper			.05 Taper		
	21mm	25mm	31mm	21mm	25mm	31mm
15	1031521	1031525	1031531	1051521	1051525	1051531



miniKUT MB2 Rotary NiTi Endodontic Negotiating Files

miniKUT MB2 Rotary NiTi Endodontic Negotiating Files (6 Packs)	
Tip Size	.05 Taper
	17mm
15	1051517

miniKUT MB2 Procedure Pack (4 Packs)	
15/.07, 15/.02 (mb2), 15/.03, 15/.05 (17mm & 25mm)	1433525



miniKUT Rotary NiTi Endodontic Orifine Shaping Files (6Packs)

Tip Size	.07 Taper
	17mm
15	1071517



miniKUT Rotary NiTi Endodontic Retreatment Files (3 Packs)

Tip Size	.06 Taper		.08 Taper	
	21mm	25mm	21mm	25mm
20	1062021	1062025		
25			1082521	1082525



Tip Size	.03 Taper	.05 Taper
15	1600315	1600515
20	1600320	1600520
25	1600325	1600525
30	1600330	1600530
35	1600335	1600535
40	1600340	1600540
45	1600345	1600545

	.03 Taper	.05 Taper
15-40	16600360	1660560



miniKUT Paper Points

Tip Size	.03 Taper	.05 Taper
15	1670315	1670515
20	1670320	1670520
25	1670325	1670525
30	1670330	1670530
35	1670335	1670535
40	1670340	1670540
45	1670345	1670545

	.03 Taper	.05 Taper
15-40	1680360	1680560



File manuali al NiTi

Le lime manuali PlanB in Nickel Titanium offrono prestazioni eccezionali. Non troverai uno strumento manuale NiTi più affidato. Un mix perfetto di flessibilità, durata e affilatura.



4020821	PlanB NiTi Endodontic Hand Files 08/.02 21mm 6 Pack
4020825	PlanB NiTi Endodontic Hand Files 08/.02 25mm 6 Pack
4020831	PlanB NiTi Endodontic Hand Files 08/.02 31mm 6 Pack
4021021	PlanB NiTi Endodontic Hand Files 10/.02 21mm 6 Pack
4021025	PlanB NiTi Endodontic Hand Files 10/.02 25mm 6 Pack
4021031	PlanB NiTi Endodontic Hand Files 10/.02 31mm 6 Pack
4021521	PlanB NiTi Endodontic Hand Files 15/.02 21mm 6 Pack
4021525	PlanB NiTi Endodontic Hand Files 15/.02 25mm 6 Pack
4021531	PlanB NiTi Endodontic Hand Files 15/.02 31mm 6 Pack
4022021	PlanB NiTi Endodontic Hand Files 20/.02 21mm 6 Pack
4022025	PlanB NiTi Endodontic Hand Files 20/.02 25mm 6 Pack
4022031	PlanB NiTi Endodontic Hand Files 20/.02 31mm 6 Pack
4022521	PlanB NiTi Endodontic Hand Files 25/.02 21mm 6 Pack
4022525	PlanB NiTi Endodontic Hand Files 25/.02 25mm 6 Pack
4022531	PlanB NiTi Endodontic Hand Files 25/.02 31mm 6 Pack
4023021	PlanB NiTi Endodontic Hand Files 30/.02 21mm 6 Pack
4023025	PlanB NiTi Endodontic Hand Files 30/.02 25mm 6 Pack
4023031	PlanB NiTi Endodontic Hand Files 30/.02 31mm 6 Pack
4023521	PlanB NiTi Endodontic Hand Files 35/.02 21mm 6 Pack
4023525	PlanB NiTi Endodontic Hand Files 35/.02 25mm 6 Pack

4023531	PlanB NiTi Endodontic Hand Files 35/.02 31mm 6 Pack
4024021	PlanB NiTi Endodontic Hand Files 40/.02 21mm 6 Pack
4024025	PlanB NiTi Endodontic Hand Files 40/.02 25mm 6 Pack
4024031	PlanB NiTi Endodontic Hand Files 40/.02 31mm 6 Pack
4024521	PlanB NiTi Endodontic Hand Files 45/.02 21mm 6 Pack
4024525	PlanB NiTi Endodontic Hand Files 45/.02 25mm 6 Pack
4024531	PlanB NiTi Endodontic Hand Files 45/.02 31mm 6 Pack
4025021	PlanB NiTi Endodontic Hand Files 50/.02 21mm 6 Pack
4025025	PlanB NiTi Endodontic Hand Files 50/.02 25mm 6 Pack
4025031	PlanB NiTi Endodontic Hand Files 50/.02 31mm 6 Pack
4025521	PlanB NiTi Endodontic Hand Files 55/.02 21mm 6 Pack
4025525	PlanB NiTi Endodontic Hand Files 55/.02 25mm 6 Pack
4025531	PlanB NiTi Endodontic Hand Files 55/.02 31mm 6 Pack
4026021	PlanB NiTi Endodontic Hand Files 60/.02 21mm 6 Pack
4026025	PlanB NiTi Endodontic Hand Files 60/.02 25mm 6 Pack
4026031	PlanB NiTi Endodontic Hand Files 60/.02 31mm 6 Pack

File manuali in acciaio inox

Non si può definire un'azienda Endo senza una linea di affidabili File manuali in acciaio inossidabile. PlanB li ha.



3020821	PlanB Stainless Steel Endodontic Hand Files 08/.02 21mm 6 Pack
3020825	PlanB Stainless Steel Endodontic Hand Files 08/.02 25mm 6 Pack
3020831	PlanB Stainless Steel Endodontic Hand Files 08/.02 31mm 6 Pack
3021021	PlanB Stainless Steel Endodontic Hand Files 10/.02 21mm 6 Pack
3021025	PlanB Stainless Steel Endodontic Hand Files 10/.02 25mm 6 Pack
3021031	PlanB Stainless Steel Endodontic Hand Files 10/.02 31mm 6 Pack
3021521	PlanB Stainless Steel Endodontic Hand Files 15/.02 21mm 6 Pack
3021525	PlanB Stainless Steel Endodontic Hand Files 15/.02 25mm 6 Pack
3021531	PlanB Stainless Steel Endodontic Hand Files 15/.02 31mm 6 Pack
3022021	PlanB Stainless Steel Endodontic Hand Files 20/.02 21mm 6 Pack
3022025	PlanB Stainless Steel Endodontic Hand Files 20/.02 25mm 6 Pack
3022031	PlanB Stainless Steel Endodontic Hand Files 20/.02 31mm 6 Pack
3022521	PlanB Stainless Steel Endodontic Hand Files 25/.02 21mm 6 Pack
3022525	PlanB Stainless Steel Endodontic Hand Files 25/.02 25mm 6 Pack
3022531	PlanB Stainless Steel Endodontic Hand Files 25/.02 31mm 6 Pack
3023021	PlanB Stainless Steel Endodontic Hand Files 30/.02 21mm 6 Pack
3023025	PlanB Stainless Steel Endodontic Hand Files 30/.02 25mm 6 Pack
3023031	PlanB Stainless Steel Endodontic Hand Files 30/.02 31mm 6 Pack
3023521	PlanB Stainless Steel Endodontic Hand Files 35/.02 21mm 6 Pack
3023525	PlanB Stainless Steel Endodontic Hand Files 35/.02 25mm 6 Pack

3023531	PlanB Stainless Steel Endodontic Hand Files 35/.02 31mm 6 Pack
3024021	PlanB Stainless Steel Endodontic Hand Files 40/.02 21mm 6 Pack
3024025	PlanB Stainless Steel Endodontic Hand Files 40/.02 25mm 6 Pack
3024031	PlanB Stainless Steel Endodontic Hand Files 40/.02 31mm 6 Pack
3024521	PlanB Stainless Steel Endodontic Hand Files 45/.02 21mm 6 Pack
3024525	PlanB Stainless Steel Endodontic Hand Files 45/.02 25mm 6 Pack
3024531	PlanB Stainless Steel Endodontic Hand Files 45/.02 31mm 6 Pack
3025021	PlanB Stainless Steel Endodontic Hand Files 50/.02 21mm 6 Pack
3025025	PlanB Stainless Steel Endodontic Hand Files 50/.02 25mm 6 Pack
3025031	PlanB Stainless Steel Endodontic Hand Files 50/.02 31mm 6 Pack
3025521	PlanB Stainless Steel Endodontic Hand Files 55/.02 21mm 6 Pack
3025525	PlanB Stainless Steel Endodontic Hand Files 55/.02 25mm 6 Pack
3025531	PlanB Stainless Steel Endodontic Hand Files 55/.02 31mm 6 Pack
3026021	PlanB Stainless Steel Endodontic Hand Files 60/.02 21mm 6 Pack
3026025	PlanB Stainless Steel Endodontic Hand Files 60/.02 25mm 6 Pack
3026031	PlanB Stainless Steel Endodontic Hand Files 60/.02 31mm 6 Pack

PlanB Endodontic Cordless Motor

Il motore PlanB (7500000) offre funzionalità complete, preimpostazioni miniKUT, una testa super piccola, ricarica wireless, un'ergonomia eccellente e un localizzatore d'apice integrato. Le funzioni di avvio e arresto automatico offrono un ulteriore livello di praticità. Ha una velocità massima di 2.500 giri/min

Una panoramica delle caratteristiche:

- Manipolo senza fili
- Localizzatore apicale integrato
- Angolo reciprocante regolabile - (intervallo di 10°; campo di regolazione: 20°- 340°)
- Contrangolo 6:1 ruotabile a 360°
- Velocità massima di 2500 giri/min
- Testina super piccola adatta al campo visivo del microscopio
- Funzioni di avvio e arresto automatici
- impostazioni preprogrammate per miniKUT e altri 10 sistemi di file in commercio personalizzabili
- Batteria al litio ricaricabile da 2000 mAh
- Carica a induzione.



L'aggiornamento atteso da tempo per le frese d'accesso.
Ingenierizzate per preparare un accesso MIE con velocità ed efficienza.

Diamantate per l'Accesso

Le diamantate per l'accesso PlanB sostituiscono ampiamente la fresa LA Axxess™*, poiché offrono l'accesso MIE ideale per la serie miniKUT e per altri marchi.

*LA Axxess è un marchio di Kerr



MIE Kit Frese

Caratteristiche del KIT FRESE MIE 3 per tipo. (9000000)

Diamantata Xtra Coarse	(9000005)
Diamantata Medium	(9000006)
Fresa arrotondata a fessura con taglio incrociato	(9000008)
Carburio N. 2	(9000004)*

*Non inclusa nel kit

