

# Saliva Ejectors & Aspirator Tubes

By **ORSING**

Orsing is one of the world's leading producers of saliva ejectors, aspirator tubes and accessories. Founded in Sweden in the 1950s, Orsing has distinguished itself for user-friendly disposable products developed with the patient's highest comfort in mind such as the Hygoformic, the world's best known mouldable saliva ejector.

**ORSING**  
directadentalgroup.com

Orsing is one of the world's leading producers of saliva ejectors, aspirator tubes and accessories. Founded in Sweden in the 1950s, Orsing has distinguished itself for user-friendly disposable products developed with the patient's highest comfort in mind such as the Hygoformic, the world's best known mouldable saliva ejector.



**J.H. ORSING AB**  
Torbornavägen 26, Helsingborg  
P.O. Box 16077, S-250 16 RÅÅ, Sweden  
Int.+46 42 29 55 00, Fax. +46 42 29 55 05  
www.orsing.com  
orsing@orsing.se

part of:  
www.directadentalgroup.com

**ORSING**  
directadentalgroup.com

## 4 reasons to use Hygoformic®

Hygoformic®

### #1

#### Adjustable

The Hygoformic is seamlessly adjustable in shape and size to any patients cheek and also works very well when treatment is needed at the back of the jaw.

### #2

#### Environmental friendly

The material in Hygoformic is a non-toxic and non-pollutive mixture of PE and PP. When the plastic is totally burnt all that will remain is water and carbon-dioxide. No phthalates – no latex.

### #3

#### High capacity and tissue friendly

The Hygoformic has a high suction capacity due to the 5 holes. The holes are placed inside the coil design to avoid contact with the tongue and soft tissues and will not become blocked.

### #4

#### Tongue holder

The Hygoformic is a combined saliva ejector and tongue holder. The coil design of the tongue holder makes Hygoformic patient friendly.

### Hygoformic®

The Hygoformic is the original mouldable saliva ejector with tongue holder. Used by dentists all over the world for more than 50 years. Produced in non-toxic and non-pollutive polyethylene and polypropylene and flanked by a complete line of adaptors to meet all the requirements of modern dental clinics.



“The world’s  
best known  
mouldable  
saliva ejector”



**Hygoformic U**  
The original mouldable saliva ejector  
Colour: white

Art.nr: U2000    Carton of 2000 pcs  
Art.nr: U2100    Carton of 20 x 100 pcs  
Art.nr: U1000    Carton of 1000 pcs

- The original mouldable saliva ejector
- With adjustable tongue holder for patient’s comfort
- Manufactured in an environmental friendly mixture of PE/PP
- Wide selection of accessories available
- Transport-economic packaging thanks to improved material and maximum utilization of space
- Certified class IIa product according to 93/42/EEC directive on medical devices

“With  
children  
in  
mind”



**Hygoformic L**  
Small saliva ejector for children  
Colour: white

Art.nr: L1000    Carton of 1000 pcs

- Smaller size saliva ejector suitable for children
- With adjustable tongue holder for patient’s comfort
- Manufactured in an environmental friendly mixture of PE/PP
- Wide selection of accessories available
- Transport-economic packaging thanks to improved material and maximum utilization of space
- Certified class IIa product according to 93/42/EEC directive on medical devices



**Soft adaptor  
for Hygoformic**

Art.no. S100  
6,5 mm  
Manufactured in PE  
Packet of 100 pcs



**Double adaptor  
for Hygoformic and other saliva ejectors**

Art.no. DA10  
Manufactured in PE  
Packet of 10 pcs



**Autoclavable hard adaptor  
for Hygoformic**

Art.no. H6 6,5 mm  
Manufactured in PP  
Packet of 10 pcs



**Autoclavable hard adaptor  
for Hygoformic**

Art.no. H6,5 6,5 mm  
Manufactured in POM  
Packet of 10 pcs



**Autoclavable hard adaptor  
for Hygoformic and other saliva ejectors**

Art.no. H11 11 mm  
Manufactured in PP  
Packet of 10 pcs




Insert an Hygoformic saliva ejector into standard ejection system. Use a suitable adaptor for Hygoformic if necessary.



The Hygoformic saliva ejector should stay in place without hurting. If the patient experiences discomfort, the Hygoformic saliva ejector should be loosened.



Stretch the Hygoformic saliva ejector to suit the patient's anatomy. The saliva ejector has a memory and can be stretched and bent to suitable shape.



Choose an aspirator tube from the wide range produced by Orsing for optimal suction and to keep the surgery area clean and dry. All Orsing's saliva ejectors and aspirator tubes can be used with a flexible AdaptaPal (Directa) for better access from all angles and minimum strain on operator's shoulders and arms.

## Disposable aspirator tubes

Orsing's wide range of disposable aspirator tubes fit any evacuation system. Besides the standard selection of models and colours, Orsing manufactures aspirator tubes in any length, design or colour to meet customers' needs.

# 4 reasons to use aspirator tubes

## #1

### Dual tips: S-shaped and 45°

The Hygovac range of aspirator tubes consists of Hygovac and the ventilated Hygovac Vent. The dual tips give the freedom of choice also in the autoclavable versions Hygovac Green and Hygovac Vent Green.

## #2

### Patient and environmentally friendly

All aspirator tubes are manufactured in an environmentally friendly PE or PP. When the plastic is totally burnt all that will remain is water and carbon-dioxide. No phthalates - no latex.



## #3

### Multi colours

The Scantube range of aspirator tubes is available in 11 colours.

## #4

### Ventilation hole

The Hygovac Vent is equipped with double holes on each end which strongly reduce the effect of getting stuck. The Scantube Vent is equipped with a rectangular ventilation hole.

### Scantube®

Plain disposable aspirator tubes with 45° + 45° cuts.  
 Length: 135 mm.  
 Manufactured in PE  
 Carton of 10 x 100 pcs

#### Art.no.


AT1000	White
BL1000	Blue
LC1000	Lilac
MG1000	Mintgreen
PU1000	Purple
PK1000	Pink
LB1000	Lightblue
OR1000	Orange
BK1000	Black
LG1000	Limegreen
YW1000	Yellow

### Angles for aspirator tube

Disposable angle for aspirator tubes.  
 Art.no. A100 Opaque  
 Manufactured in PE.  
 Packet of 100 pcs

#### Art.no.

A100	Opaque
A100 x 50	Opaque
A100LB	Lightblue
A100MG	Mintgreen
A100OR	Orange
A100PU	Purple



45°

S-shaped

### Hygovac®

Disposable aspirator tube with one traditional 45° end and one S-shaped end for increased patient comfort.  
 Colour: White. Length: 140 mm.  
 Manufactured in PE.  
 Art.no. V2800  
 Carton of 28 x 100 pcs  
 Art.no. V1000  
 Carton of 10 x 100 pcs



45°

S-shaped

### Hygovac® Vent

Double-ventilated disposable aspirator tube with one traditional 45° end and one S-shaped end. The double holes on each end strongly reduce the effect of getting stuck.  
 Colour: White with blue stripes.  
 Length: 140 mm.  
 Manufactured in PE.  
 Art.no. V2800V  
 Carton of 28 x 100 pcs  
 Art.no. V1000V  
 Carton of 10 x 100 pcs



45°

ventilation hole

S-shaped

### Scantube® Vent

Disposable aspirator tube with S-shaped ends and a rectangular ventilation hole. Length: 145 mm.  
 Manufactured in PE.  
 Art.no. AV1000 White  
 AV1000LB Lightblue  
 AV1000MG Mintgreen  
 AV1000OR Orange  
 AV1000PU Purple  
 AV1000BL Blue  
 Carton of 10 x 100 pcs



### Angles for aspirator tube

Disposable angle for aspirator tubes. Manufactured in PE.  
 Art.no. A100 Opaque  
 Packet of 100 pcs  
 Other colours on request.

## Autoclavable aspirator tubes & accessories

Orsing's line of autoclavable products. Both the long and the short autoclave program are suitable (max 140° C). The tips Hygosurge, Hygojet and Hygotip as well as the disposable angle fit all our aspirator tubes.







## 4 reasons to use Hygofix®

### #1

#### One piece design

The Hygofix has an integrated tip.  
No loose tip that can fall off.

### #2

#### Pre-bent

The Hygofix is prebent in an optimal fixed design for total patient comfort.



### #3

#### High capacity and tissue friendly

The Hygofix has a high suction capacity and is tissue friendly. The design of the integrated tip ensures patient comfort and non clogging operation.

### #4

#### Environmental friendly and manufactured in one material

The Hygofix is manufactured in one material only a non-toxic and non-pollutive PP. No phthalates – no latex.

#### Hygofix®

Prebent saliva ejector with integrated tip to fit 6,5 mm system.  
Manufactured in one piece and one material.

Colour: White.

Manufactured in PP.

Carton of 10x100 pcs

Art.no. HF1100



## 4 reasons to use Hygosurge® XL

### #1 Ideal for surgical treatments

The Hygosurge XL works directly in the 11 mm system and is ideal for surgical treatments.

### #2 Environmental friendly

The Hygosurge XL is manufactured in a non-toxic and non-pollutive PP. No phthalates – No latex



### #3 Autoclavable

The Hygosurge XL can be autoclaved on both the long and short program on max 140° C before use.

### #4 Precision tip

The Hygosurge XL has an apex of 3,6 mm and for micro oral surgery and fine endodontic work we recommend to use a combination with Microtip with an apex of 1,7 mm.

#### Hygosurge® XL

Autoclavable surgical aspirator to fit 11 mm system. Apex: 3,6 mm. Length: 165 mm. Colour: White. Autoclavable. Manufactured in PP. Art.no. HS25XL. Packet of 25 pcs



**Hygovac® Green**

Autoclavable aspirator tube with one traditional 45° end and one S-shaped end for increased patient comfort. Autoclavable. Colour: Green. Length: 140 mm. Manufactured in PP. Art.no. V1000G Carton of 10 x 100 pcs



**Hygovac® Vent Green**

The new double-ventilated autoclavable aspirator tube with one traditional 45° end and one S-shaped end. The double holes on each end strongly reduce the effect of getting stuck. Colour: White with green stripes. Length: 140 mm. Manufactured in PP. Autoclavable. Art.no. V25VG Packet of 25 pcs



**Hygosurge®**

This surgical tip has an apex of only 3,6 mm, which naturally provides increased possibilities. It is perfect to use when dealing with narrow operations, root fillings and implants. Colour: White. Autoclavable. Manufactured in PP. Art.no. HS25 Packet of 25 pcs



**Hygojet®**

The Hygojet aspirator tip is silent to use and easy to handle. Its design enables the dentist to reach by holding back the tongue as well as the cheek. Colour: White. Autoclavable. Manufactured in PP. Art.no. HJ25 Packet of 25 pcs



**Hygotip®**

The Hygotip is comfortable and protective to use and thanks to its smallness, very suitable when treating children. Colour: White. Autoclavable. Manufactured in PP. Art.no. TP25 Packet of 25 pcs



**Hygotip® XL**

Autoclavable aspirator to fit 16 mm system. Length: 145 mm. Colour: White. Manufactured in PP. Packet of 10 pcs Art.no. TP10XL



**Microtip®**

A tiny blue autoclavable tip that turns the Hygosurge and Hygosurge XL into a unique suction device for micro oral surgery and fine endodontic work. Apex: 1,7 mm. Colour: Lightblue. Autoclavable. Manufactured in PP. Art.no. MT25 Packet of 25 pcs



**H16**

Adaptor for Hygovac, other aspirator tubes and Hygosurge XL to fit 16 mm system. Colour: White. Autoclavable. Manufactured in PP. Art.no. H11-16 Packet of 10 pcs



**H16XL**

Adaptor for Hygovac, other aspirator tubes and Hygosurge XL to fit CATTANI 16 mm system. Colour: White. Autoclavable. Manufactured in POM. Art.no. H16XL Packet of 10 pcs



**H11-16**

Reduction adaptor to fit 11 mm system for Hygotip XL and other aspirator tubes for 16 mm system. Colour: White. Autoclavable. Manufactured in POM. Art.no. H11-16 Packet of 10 pcs



**Safety container**

For disposing of hazardous waste. 2200 ml. Manufactured in PP. Art.no. SC25 Carton of 25 pcs

# Product Numbers

## Saliva Ejectors Hygoformic® U

Universal saliva ejector. Colour: White. Manufactured in PE/PP.

Prod No	Product Description	Quantity
U2000	Saliva Ejectors Hygoformic U	Price per carton of 1 x 2000 pcs. 60 adaptors included. 40 cartons = full pallet
U2100	Saliva Ejectors Hygoformic U	Price per carton of 20 x 100 pcs. 5 adaptors per bag. 40 cartons = full pallet
U1000	Saliva Ejectors Hygoformic U	Price per carton of 1 x 1000 pcs. 30 adaptors included. 80 cartons = full pallet

## Saliva Ejectors Hygoformic® L

Small for children. Colour: White. Manufactured in PE/PP.

Prod No	Product Description	Quantity
L1000	Saliva Ejectors Hygoformic L	Price per carton of 1 x 1000 pcs. Price per carton of 1 x 1000 pcs. 30 adaptors. 80 cartons = full pallet included

## Hygofix®

Prebended saliva ejector with integrated tip to fit 6,5 mm system. Manufactured in one piece and one material. Colour: White. Manufactured in PP.

Prod No	Product Description	Quantity
HF1100	Hygofix	Price per carton of 10 x 100 pcs.

## Adaptors for Hygoformic® and other saliva ejectors

To fit into standard ejection systems.

Prod No	Product Description	Quantity
S100	Soft adaptor, 6,5 mm PE	Per packet of 100 pcs
S100 x 100	Carton with 100 packets	Per packet of 100 pcs
H6	Hard adaptor, 6,5 mm PP	Per packet of 10 pcs
H6,5	Hard adaptor, 6,5 mm POM	Per packet of 10 pcs
H11	Hard adaptor, 11 mm PP	Per packet of 10 pcs
DA10	Double adaptor PE	Per packet of 10 pcs
T10	T-part, double use PE/PP	Per packet of 10 pcs

## Aspirator tubes Hygovac®

Aspirator tube with S-cut + 45° cut. Colour: White. Manufactured in PE.

Prod No	Product Description	Quantity
V2800	Hygovac	Price per carton of 28 x 100 pcs. 32 cartons = full pallet
V1000	Hygovac	Price per carton of 10 x 100 pcs. 80 cartons = full pallet

Ventilated aspirator tube with S-cut + 45° cut. Colour: White with blue stripes. Manufactured in PE.

Prod No	Product Description	Quantity
V2800V	Hygovac Vent	Price per carton of 28 x 100 pcs. 32 cartons = full pallet
V1000V	Hygovac Vent	Price per carton of 10 x 100 pcs. 80 cartons = full pallet
V1000G	Hygovac Green, Colour: Green	Price per carton of 10 x 100 pcs. 80 cartons = full pallet
V25VG	Hygovac Vent Green, white with green stripes	Price per packet of 25 pcs.

## Scantube®

Aspirator tubes with 45° + 45° cuts. Length: 135 mm.

Prod No	Product Description	Quantity
AT1000	White	Price per carton of 10 x 100 pcs
LG1000	Limegreen	80 cartons = full pallet
YW1000	Yellow	Price per carton of 10 x 100 pcs
LC1000	Lilac	80 cartons = full pallet
BK1000	Black	Price per carton of 10 x 100 pcs.
PK1000	Pink	80 cartons = full pallet
LB1000	Light-blue	Price per carton of 10 x 100 pcs.
MG1000	Mintgreen	80 cartons = full pallet
OR1000	Orange	Price per carton of 10 x 100 pcs.
PU1000	Purple	80 cartons = full pallet
BL1000	Blue	Price per carton of 10 x 100 pcs.

## Accessories

For Hygovac, Hygovac Vent and our other aspirator tubes

Prod No	Product Description	Quantity
HS25	Hygosurge, autoclavable aspirator tips. Colour: White.	Price per packet of 25 pcs
HS25 x 100		Price per carton of 100 packets
HJ25	Hygojet Autoclavable aspirator tips. Colour: White.	Price per packet of 25 pcs
HJ25 x 100		Price per carton of 100 packets
TP25	Hygotip, Autoclavable aspirator tips. Colour: White.	Price per packet of 25 pcs
TP25 x 100		Price per carton of 100 packets
MT25	Microtip, colour: lightblue.	Price per packet of 25 pcs
H16	Adaptor for Hygovac, Hygovac Vent and Hygosurge XL	Price per packet of 10 pcs
H16XL	Adaptor for Hygovac, Hygovac Vent and Hygosurge XL to fit CATTANI 16 mm system.	Price per packet of 10 pcs
H11-16	Reduction Adaptor	Price per packet of 10 pcs
A100	Disposable angle for aspirator tube Opaque	Price per packet of 100 pcs
A100 x 50	Disposable angle for aspirator tube Opaque	Price per carton of 50 x 100 pcs
A100LB	Disposable angle for aspirator tube Lightblue	Price per packet of 100 pcs
A100MG	Disposable angle for aspirator tube Mintgreen	Price per packet of 100 pcs
A100OR	Disposable angle for aspirator tube Orange	Price per packet of 100 pcs
A100PU	Disposable angle for aspirator tube Purple	Price per packet of 100 pcs
HS25XL	Hygosurge XL surgical aspirator Autoclavable aspirator to fit 11 mm system. Colour: White. Apex 3,6 mm. Length 165 mm.	Price per packet of 25 pcs
TP10XL	Hygotip XL Autoclavable aspirator to fit 16 mm system. Colour: White. Length 150 mm.	Price per packet of 10 pcs.
HE25ST	Hygoscreen Straight Autoclavable aspirator to fit 11 mm system. Colour: White. Length 155 mm.	Price per packet of 25 pcs.
PT100	Plastic Tubing, Colour: Transparent.	Price per roll of 100 m.
SC25	Safety Container, 2,2 litres. Colour: Yellow	Price per carton of 25 pcs with lids.w

